

Stem Cell Populations Detected by HemoGenix® Assays

Clonal, Methylcellulose
Colony-Forming Cell
Differentiation Assays
Available




- ColonyGro™:**
Traditional CFU assay
- CAMEO™-4:**
Miniaturized CFU assay
- CAMEO™-96:**
Standardized CFU assay that detects proliferation and differentiation

Instrument-Based,
Non-Methylcellulose
Stem Cell Proliferation
Assays Available



- HALO®-96 Research:**
ATP bioluminescence readout
- HemoFLUOR™-96 Research:**
Fluorescence readout
- HemoLIGHT™-96 Research:**
Absorbance readout



All clonal and instrument-based assays are available for:

- Human
- Non-human primate
- Horse
- Pig
- Sheep
- Dog
- Rat
- Mouse

Serum-free, instrument-based assays are available for:

- Human
- Non-human primate
- Mouse